

scikit-surgerybk Documentation

Thomas Dowrick

Feb 24, 2022

Contents

1 Developing	3
2 Installing	5
3 Licensing and copyright	7
4 Acknowledgements	9
Python Module Index	15
Index	17



Author: Thomas Dowrick

scikit-surgerybk is a Python library for interfacing with BK ultrasound machines, and is part of the [SciKit-Surgery](#) software project, developed at the [Wellcome EPSRC Centre for Interventional and Surgical Sciences](#), part of [University College London \(UCL\)](#).

scikit-surgerybk supports Python 2.7 and Python 3.6/3.7/3.8.

CHAPTER 1

Developing

1.1 Cloning

You can clone the repository using the following command:

```
git clone https://github.com/SciKit-Surgery/scikit-surgerybk
```

1.2 Running tests

Pytest is used for running unit tests:

```
pip install pytest
python -m pytest
```

1.3 Linting

This code conforms to the PEP8 standard. Pylint can be used to analyse the code:

```
pip install pylint
pylint --rcfile=tests/pylintrc sksurgerybk
```


CHAPTER 2

Installing

You can pip install directly from the repository as follows:

```
pip install git+https://github.com/SciKit-Surgery/scikit-surgerybk
```

2.1 Contributing

Please see the contributing guidelines.

2.2 Useful links

- Source code repository
- Documentation

CHAPTER 3

Licensing and copyright

Copyright 2019 University College London. scikit-surgerybk is released under the BSD-3 license. Please see the [license file](#) for details.

CHAPTER 4

Acknowledgements

Supported by Wellcome and EPSRC.

4.1 Requirements for scikit-surgerybk

This is the software requirements file for scikit-surgerybk, part of the SNAPPY project. The requirements listed below should define what scikit-surgerybk does. Each requirement can be matched to a unit test that checks whether the requirement is met.

4.1.1 Requirements

ID	Description	Test
0000	Module has a help page	pylint, see tests/pylint.rc and tox.ini
0001	Functions are documented	pylint, see tests/pylint.rc and tox.ini
0002	Package has a version number	No test yet, handled by git.

4.2 latest

4.2.1 sksurgerybk package

Subpackages

sksurgerybk.interface package

Submodules

sksurgerybk.interface.bk5000 module

This module sets the connection to the BK scanner

class `sksurgerybk.interface.bk5000.BK5000(timeout, frames_per_second)`

Bases: `object`

This class sets the TCP connection with the BK scanner

clear_bytes_in_buffer(start, end)

Clear a set of bytes in bytearray buffer

Parameters

- **start** (`integer`) – Start index
- **end** (`integer`) – End integer

connect_to_host(address, port)

Connects the client to the host/serverself.

Implements a try/except block to catch potential errors.

Parameters: address(`string`): the IP address port(`integer`): the port

decode_image()

Process the stream of data received from the BK5000 and convert it into a numpy array which represents the ultrasound image.

Control bytes are 1, 4 and 27. Flipped control bytes (1s complement of control bytes) are 254, 251, 228. Any time a flipped control bytes occurs after a 27, the value should be flipped and the preceding 27 deleted. See page 9 of 142 in BK doc PS12640-44 for further details.

disconnect_from_host()

Disconnects the client from the host.

If the socket is already closed, a `recv()` call will throw an error. If it doesn't, we can close the socket.

enable_rgb_output()

The 'QUERY:GRAB_FRAME "ON"' gets the BK to stream greyscale data (e.g. 640 x 480 x 1),

some applications might want this in RGB format (640 x 480 x 3)

find_first_a_not_preceded_by_b(start_pos, a, b)

Find the first instance of 'a' in an array that isn't preceded by 'b'

Parameters

- **start_pos** (`integer`) – Index in array to begin search at
- **a** (`integer`) – Value to find
- **b** (`integer`) – Value not to precede a

Returns Index of first a not preceded by b, -1 if none found

Return type `integer`

generate_command_message(message)

Append 0x01 and 0x04 to start/end of message before sending

Parameters: message(`string`): the message to be sent

get_frame()

Get the next frame from the BK5000.

parse_scanarea_message()

Separate the scanarea response message into the separate components.

Message has format: <StartLineX>,<StartLineY>,<StartLineAngle>,<StartDepth>, <StopLineX>,<StopLineY>,<StopLineAngle>,<StopDepth>

Example message: DATA:B_GEOMETRY_SCANAREA:A 0.0017218,-0.000171398, 1.37236,0,-0.00174855,-0.000176821,1.77236,0.0203479;

parse_win_size_message(*message*)

Extract the size of the US window from the response message

Message has format “DATA:US_WIN_SIZE 640,480;”

Parameters: message(string): the received message

query_scanarea()

Query the BK5000 for the scan area. (width/height/scanning depth etc.)

query_win_size()

Query the BK5000 for the window/image size

receive_image()

Scan the incoming data stream to find the start and end of the image data.

See BK doc PS12640-44 for further details.

receive_response_message(*expected_size=1024*)

Receive a message

Stores it under the data class member

Parameters: expected_size(int): the receive message size in bytes

request_stop()

Set the appropriate class member

send_command_message(*message*)

Send a message through the socket.

Implements a couple of checks to verify the message has been sent correctly.

Parameters: message(string): the message to be sent

start_streaming()

Send a message to start the streaming

stop_streaming()

Send a message to stop the streaming. send_command_message and receive_response_message will throw errors if there is a problem with the socket connection.

```
class sksurgerybk.interface.bk5000.BKOpenCV (TCP_IP='128.16.0.3', TCP_PORT=7915,  
 TIMEOUT=5, FPS=25)
```

Bases: object

Display BK data using OpenCV.

start()

Start acquisition/streaming.

```
class sksurgerybk.interface.bk5000.BKpyIGTLink (TCP_IP='128.16.0.3',  
 TCP_PORT=7915, TIMEOUT=5,  
 FPS=8)
```

Bases: object

Send BK data over OpenIGTLink.

```
start()
    Start acquisition/streaming.

stop()
    Stop acquisition/streaming.

sksurgerybk.interface.bk5000.main()
    Entry point for OpenCV/pyIGTLink.
```

Module contents

sksurgerybk.pyigtlink package

Submodules

sksurgerybk.pyigtlink.pyIGTLink module

Created on Tue Nov 3 19:17:05 2015

@author: Daniel Hoyer Iversen

```
class sksurgerybk.pyigtlink.pyIGTLink.ImageMessage(image, spacing=[1, 1, 1], timestamp=None, device_name="")
    Bases: sksurgerybk.pyigtlink.pyIGTLink.MessageBase

    pack_body()

class sksurgerybk.pyigtlink.pyIGTLink.ImageMessageMatlab(image, dim, spacing=[1, 1, 1], timestamp=None)
    Bases: sksurgerybk.pyigtlink.pyIGTLink.ImageMessage

class sksurgerybk.pyigtlink.pyIGTLink.MessageBase
    Bases: object

    message

    get_binary_body()
    get_binary_message()
    get_body_pack_size()
    is_valid()

    pack()
    pack_body()

class sksurgerybk.pyigtlink.pyIGTLink.PyIGTLink(port=18905, localServer=True, iface='eth0')
    Bases: socketserver.TCPServer

    For streaming data over TCP with IGTLINK

    add_message_to_send_queue(message, wait=False)
        Returns True if sucessfull

    close_server()
        Will close connection and shutdown server

    get_ip_adress()
    get_port_no()
```

```
is_connected()
update_connected_status(val)

class sksurgerybk.pyigtlink.pyIGTLink.TCPRequestHandler(request, client_address,
server)
    Bases: socketserver.BaseRequestHandler

    Help class for PyIGTLink

    handle()

class sksurgerybk.pyigtlink.pyIGTLink.TransformMessage(tform, timestamp=None, device_name=‘’)
    Bases: sksurgerybk.pyigtlink.pyIGTLink.MessageBase

    pack_body()
```

Module contents

Module contents

scikit-surgerybk

4.2.2 work_in_progress module

- modindex
- genindex
- search

Python Module Index

i

`sksurgerybk.interface`, 12
`sksurgerybk.interface.bk5000`, 10

p

`sksurgerybk.pyigtlink`, 13
`sksurgerybk.pyigtlink.pyIGTLink`, 12

s

`sksurgerybk`, 13

Index

A

add_message_to_send_queue()
(*sksurgerybk.pyiglink.pyIGTLink.PyIGTLink method*), 12

B

BK5000 (*class in sksurgerybk.interface.bk5000*), 10
BKOpenCV (*class in sksurgerybk.interface.bk5000*), 11
BKpyIGTLink (*class in sksurgerybk.interface.bk5000*), 11

C

clear_bytes_in_buffer()
(*sksurgerybk.interface.bk5000.BK5000 method*), 10
close_server() (*sksurgerybk.pyiglink.pyIGTLink.PyIGTLink method*), 12
connect_to_host()
(*sksurgerybk.interface.bk5000.BK5000 method*), 10

D

decode_image() (*sksurgerybk.interface.bk5000.BK5000 method*), 10
disconnect_from_host()
(*sksurgerybk.interface.bk5000.BK5000 method*), 10

E

enable_rgb_output()
(*sksurgerybk.interface.bk5000.BK5000 method*), 10

F

find_first_a_not_preceded_by_b()
(*sksurgerybk.interface.bk5000.BK5000 method*), 10

G

generate_command_message()
(*sksurgerybk.interface.bk5000.BK5000 method*), 10
get_binary_body()
(*sksurgerybk.pyiglink.pyIGTLink.MessageBase method*), 12
get_binary_message()
(*sksurgerybk.pyiglink.pyIGTLink.MessageBase method*), 12
get_body_pack_size()
(*sksurgerybk.pyiglink.pyIGTLink.MessageBase method*), 12
get_frame() (*sksurgerybk.interface.bk5000.BK5000 method*), 10
get_ip_adress() (*sksurgerybk.pyiglink.pyIGTLink.PyIGTLink method*), 12
get_port_no() (*sksurgerybk.pyiglink.pyIGTLink.PyIGTLink method*), 12

H

handle() (*sksurgerybk.pyiglink.pyIGTLink.TCPRequestHandler method*), 13

|

ImageMessage (*class in sksurgerybk.pyiglink.pyIGTLink*), 12
ImageMessageMatlab (*class in sksurgerybk.pyiglink.pyIGTLink*), 12
is_connected() (*sksurgerybk.pyiglink.pyIGTLink.PyIGTLink method*), 13
is_valid() (*sksurgerybk.pyiglink.pyIGTLink.MessageBase method*), 12

M

main() (*in module sksurgerybk.interface.bk5000*), 12
MessageBase (*class in sksurgerybk.pyiglink.pyIGTLink*), 12

P

pack () (*sksurgerybk.pyigtlink.pyIGTLink.MessageBase* method), 12
pack_body () (*sksurgerybk.pyigtlink.pyIGTLink.ImageMessage* method), 12
pack_body () (*sksurgerybk.pyigtlink.pyIGTLink.MessageBase* method), 12
pack_body () (*sksurgerybk.pyigtlink.pyIGTLink.TransformMessage* method), 13
parse_scanarea_message () (*sksurgerybk.interface.bk5000.BK5000* method), 10
parse_win_size_message () (*sksurgerybk.interface.bk5000.BK5000* method), 11
PyIGTLink (*class in sksurgerybk.pyigtlink.pyIGTLink*), 12

T

TCPRequestHandler (*class in sksurgerybk.pyigtlink.pyIGTLink*), 13
MessagesformMessage (*class in sksurgerybk.pyigtlink.pyIGTLink*), 13
U
connected_status () (*sksurgerybk.pyigtlink.PyIGTLink* method), 13

Q

query_scanarea () (*sksurgerybk.interface.bk5000.BK5000* method), 11
query_win_size () (*sksurgerybk.interface.bk5000.BK5000* method), 11

R

receive_image () (*sksurgerybk.interface.bk5000.BK5000* method), 11
receive_response_message () (*sksurgerybk.interface.bk5000.BK5000* method), 11
request_stop () (*sksurgerybk.interface.bk5000.BK5000* method), 11

S

send_command_message () (*sksurgerybk.interface.bk5000.BK5000* method), 11
sksurgerybk (*module*), 13
sksurgerybk.interface (*module*), 12
sksurgerybk.interface.bk5000 (*module*), 10
sksurgerybk.pyigtlink (*module*), 13
sksurgerybk.pyigtlink.pyIGTLink (*module*), 12
start () (*sksurgerybk.interface.bk5000.BKOpenCV* method), 11
start () (*sksurgerybk.interface.bk5000.BKpyIGTLink* method), 11
start_streaming () (*sksurgerybk.interface.bk5000.BK5000* method), 11
stop () (*sksurgerybk.interface.bk5000.BKpyIGTLink* method), 12
stop_streaming () (*sksurgerybk.interface.bk5000.BK5000* method), 11